

For Research Use Only

KDM3A, JMJD1A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83571-4-PBS

Basic Information

Catalog Number: 83571-4-PBS	GenBank Accession Number: BC038297	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55818	CloneNo.: 240532D2
Source: Rabbit	UNIPROT ID: Q9Y4C1	
Isotype: IgG	Full Name: lysine (K)-specific demethylase 3A	
Immunogen Catalog Number: AG28687	Calculated MW: 150 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00565-1, KDM3A, JMJD1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83571-2-PBS. Detection antibody: 83571-4-PBS. Standard: Ag28687. Range: 1.56-100 ng/mL.