

For Research Use Only

Mouse VEGFA164 Recombinant antibody, PBS Only (Detector)



Catalog Number:83572-1-PBS

Basic Information

Catalog Number:

83572-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001287057

GeneID (NCBI):

22339

UNIPROT ID:

Q00731-2

Full Name:

vascular endothelial growth factor A

Calculated MW:

22kd

Purification Method:

Protein A purification

CloneNo.:

240415A6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00561-3, MOUSE VEGFA164 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83572-3-PBS. Detection antibody: 83572-1-PBS. Standard: Eg0277. Range: 0.625-80 ng/mL.