

For Research Use Only

TNFB Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83573-1-PBS

Basic Information

Catalog Number: 83573-1-PBS	GenBank Accession Number: BC034729	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4049	CloneNo.: 240393C2
Source: Rabbit	UNIPROT ID: P01374	
Isotype: IgG	Full Name: lymphotoxin alpha (TNF superfamily, member 1)	
	Calculated MW: 205 aa, 22 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00562-2, TNF beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83573-4-PBS. Detection antibody: 83573-1-PBS. Standard: Eg0836. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00562-3, TNF beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83573-1-PBS. Detection antibody: 83573-2-PBS. Standard: Eg0836. Range: 0.313-40 ng/mL.