

For Research Use Only

CD320 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83574-2-PBS

Basic Information

Catalog Number: 83574-2-PBS	GenBank Accession Number: BC000668	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51293	CloneNo.: 240500H5
Source: Rabbit	UNIPROT ID: Q9NPF0	
Isotype: IgG	Full Name: CD320 molecule	
Immunogen Catalog Number: AG21445	Calculated MW: 29 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00564-2, CD320 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83574-2-PBS. Detection antibody: 83574-4-PBS. Standard: Ag21445. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00564-3, CD320 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83574-3-PBS. Detection antibody: 83574-2-PBS. Standard: Ag21445. Range: 1.25-40 ng/mL