

For Research Use Only

# CD320 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83574-4-PBS

## Basic Information

Catalog Number:

83574-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21445

GenBank Accession Number:

BC000668

GeneID (NCBI):

51293

UNIPROT ID:

Q9NPF0

Full Name:

CD320 molecule

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

240500D2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00564-2, CD320 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83574-2-PBS. Detection antibody: 83574-4-PBS. Standard: Ag21445. Range: 0.625-40 ng/mL.