

For Research Use Only

S100A9 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83578-3-PBS

Basic Information

Catalog Number: 83578-3-PBS	GenBank Accession Number: BC047681	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6280	CloneNo.: 240539D5
Source: Rabbit	UNIPROT ID: P06702	
Isotype: IgG	Full Name: S100 calcium binding protein A9	
Immunogen Catalog Number: AG25764	Calculated MW: 13 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00569-1, S100A9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83578-3-PBS. Detection antibody: 83578-1-PBS. Standard: Ag25764. Range: 0.313-40 ng/mL