

For Research Use Only

# SEMA3A Recombinant antibody, PBS Only (Detector)



Catalog Number:83579-2-PBS

## Basic Information

Catalog Number:

83579-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27173

GenBank Accession Number:

BC111416

GeneID (NCBI):

10371

UNIPROT ID:

Q14563

Full Name:

sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3A

Calculated MW:

771 aa, 89 kDa

Purification Method:

Protein A purification

CloneNo.:

240509G1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00552-1, SEMA3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83579-1-PBS. Detection antibody: 83579-2-PBS. Standard: Ag27173. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00552-2, SEMA3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83579-3-PBS. Detection antibody: 83579-2-PBS. Standard: Ag27173. Range: 0.625-40 ng/mL.