

For Research Use Only

SEMA3A Recombinant antibody, PBS Only (Capture)



Catalog Number:83579-3-PBS

Basic Information

Catalog Number: 83579-3-PBS	GenBank Accession Number: BC111416	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10371	CloneNo.: 240509F11
Source: Rabbit	UNIPROT ID: Q14563	
Isotype: IgG	Full Name: sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3A	
Immunogen Catalog Number: AG27173	Calculated MW: 771 aa, 89 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00552-2, SEMA3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83579-3-PBS. Detection antibody: 83579-2-PBS. Standard: Ag27173. Range: 0.625-40 ng/mL.