

For Research Use Only

SEMA3A Recombinant antibody, PBS Only



Catalog Number: 83579-5-PBS

Basic Information

Catalog Number:

83579-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27173

GenBank Accession Number:

BC111416

GeneID (NCBI):

10371

UNIPROT ID:

Q14563

Full Name:

sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3A

Calculated MW:

771 aa, 89 kDa

Purification Method:

Protein A purification

CloneNo.:

240509E3

Applications

Tested Applications:

ELISA, WB

Species Specificity:

rat, human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83579-5-RR (SEMA3A antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83579-5-PBS in a different storage buffer formulation.