

For Research Use Only

CCDC94 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83580-3-PBS

Basic Information

Catalog Number: 83580-3-PBS	GenBank Accession Number: BC000561	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55702	CloneNo.: 240517D1
Source: Rabbit	UNIPROT ID: Q9BW85	
Isotype: IgG	Full Name: coiled-coil domain containing 94	
Immunogen Catalog Number: AG13764	Calculated MW: 323 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00553-3, CCDC94 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83580-1-PBS. Detection antibody: 83580-3-PBS. Standard: Ag13764. Range: 0.78-100 ng/mL.