

For Research Use Only

DHRS3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83581-1-PBS

Basic Information

Catalog Number:

83581-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7589

GenBank Accession Number:

BC002730

GeneID (NCBI):

9249

UNIPROT ID:

O75911

Full Name:

dehydrogenase/reductase (SDR family) member 3

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

240545C4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00554-2, DHRS3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83581-2-PBS. Detection antibody: 83581-1-PBS. Standard: Ag7589. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00554-3, DHRS3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83581-3-PBS. Detection antibody: 83581-1-PBS. Standard: Ag7589. Range: 0.313-40 ng/mL