

For Research Use Only

TRIM15 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83582-4-PBS

Basic Information

Catalog Number:

83582-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4533

GenBank Accession Number:

BC038585

GeneID (NCBI):

89870

UNIPROT ID:

Q9C019

Full Name:

tripartite motif-containing 15

Calculated MW:

465 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

240565F11

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00555-3, TRIM15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83582-2-PBS. Detection antibody: 83582-4-PBS. Standard: Ag4533. Range: 0.78-100 ng/mL