

For Research Use Only

GLIPR1L1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83584-1-PBS

Basic Information

Catalog Number: 83584-1-PBS	GenBank Accession Number: BC014603	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 256710	CloneNo.: 240594D3
Source: Rabbit	UNIPROT ID: Q6UWM5	
Isotype: IgG	Full Name: GLI pathogenesis-related 1 like 1	
Immunogen Catalog Number: AG14852	Calculated MW: 242 aa, 27 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00557-2, GLIPR1L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83584-2-PBS. Detection antibody: 83584-1-PBS. Standard: Ag14852. Range: 0.78-100 ng/mL.