

For Research Use Only

GLIPR1L1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83584-2-PBS

Basic Information

| | | |
|---|--|---|
| Catalog Number: 83584-2-PBS | GenBank Accession Number: BC014603 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 256710 | CloneNo.: 240594F4 |
| Source: Rabbit | UNIPROT ID: Q6UWM5 | |
| Isotype: IgG | Full Name: GLI pathogenesis-related 1 like 1 | |
| Immunogen Catalog Number: AG14852 | Calculated MW: 242 aa, 27 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00557-1, GLIPR1L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83584-2-PBS. Detection antibody: 83584-3-PBS. Standard: Ag14852. Range: 0.78-100 ng/mL.



Cytometric bead array standard curve of MP00557-2, GLIPR1L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83584-2-PBS. Detection antibody: 83584-1-PBS. Standard: Ag14852. Range: 0.78-100 ng/mL.