

For Research Use Only

# LEFTY1 Recombinant antibody, PBS Only



Catalog Number: 83586-1-PBS

## Basic Information

Catalog Number:

83586-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4925

GenBank Accession Number:

BC027883

GeneID (NCBI):

10637

UNIPROT ID:

O75610

Full Name:

left-right determination factor 1

Calculated MW:

366 aa, 41 kDa

Purification Method:

Protein A purification

CloneNo.:

240346E11

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83586-1-RR (LEFTY1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83586-1-PBS in a different storage buffer formulation.