

For Research Use Only

# PEBP4 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83589-3-PBS



## Basic Information

Catalog Number:

83589-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9710

GenBank Accession Number:

BC020779

GeneID (NCBI):

157310

UNIPROT ID:

Q96596

Full Name:

phosphatidylethanolamine-binding protein 4

Calculated MW:

223 aa, 25 kDa

Purification Method:

Protein A purification

CloneNo.:

240547B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00574-1, PEBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83589-1-PBS. Detection antibody: 83589-3-PBS. Standard: Ag9710. Range: 3.13-200 ng/mL



Cytometric bead array standard curve of MP00574-3, PEBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83589-3-PBS. Detection antibody: 83589-4-PBS. Standard: Ag9710. Range: 3.13-200 ng/mL