

For Research Use Only

PEBP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83589-4-PBS



Basic Information

Catalog Number: 83589-4-PBS	GenBank Accession Number: BC020779	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 157310	CloneNo.: 240547H2
Source: Rabbit	UNIPROT ID: Q96596	
Isotype: IgG	Full Name: phosphatidylethanolamine-binding protein 4	
Immunogen Catalog Number: AG9710	Calculated MW: 223 aa, 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00574-2, PEBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83589-2-PBS. Detection antibody: 83589-4-PBS. Standard: Ag9710. Range: 1.56-200 ng/mL



Cytometric bead array standard curve of MP00574-3, PEBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83589-3-PBS. Detection antibody: 83589-4-PBS. Standard: Ag9710. Range: 3.13-200 ng/mL