

For Research Use Only

INTS4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83590-2-PBS

Basic Information

Catalog Number:

83590-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9116

GenBank Accession Number:

BC009995

GeneID (NCBI):

92105

UNIPROT ID:

Q96HW7

Full Name:

integrator complex subunit 4

Calculated MW:

963aa, 108 kDa; 505aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

240579H7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00580-1, INTS4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83590-5-PBS. Detection antibody: 83590-2-PBS. Standard: Ag9116. Range: 0.78-100 ng/mL