

For Research Use Only

INTS4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83590-4-PBS

Basic Information

Catalog Number: 83590-4-PBS	GenBank Accession Number: BC009995	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 92105	CloneNo.: 240579B10
Source: Rabbit	UNIPROT ID: Q96HW7	
Isotype: IgG	Full Name: integrator complex subunit 4	
Immunogen Catalog Number: AG9116	Calculated MW: 963aa, 108 kDa; 505aa, 57 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00580-2, INTS4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83590-4-PBS. Detection antibody: 83590-3-PBS. Standard: Ag9116. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00580-3, INTS4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83590-4-PBS. Detection antibody: 83590-1-PBS. Standard: 6.25-100 ng/mL. Range: 6.25-100 ng/mL