

For Research Use Only

# Neurogranin Recombinant antibody, PBS Only (Capture)



Catalog Number: 83593-3-PBS

## Basic Information

|                                     |  |  |
|-------------------------------------|--|--|
| Catalog Number:<br>83593-3-PBS      | GenBank Accession Number:<br>BC002835                          | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                    | GeneID (NCBI):<br>4900   | CloneNo.:<br>240526F9                          |
| Source:<br>Rabbit                   | UNIPROT ID:<br>Q92686  |  |
| Isotype:<br>IgG                     | Full Name:<br>neurogranin (protein kinase C<br>substrate, RC3) |  |
| Immunogen Catalog Number:<br>AG0676 | Calculated MW:<br>7.6 kDa                                      |  |

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP00606-2, NRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83593-3-PBS. Detection antibody: 83593-2-PBS. Standard: Ag0676. Range: 0.156-20 ng/mL