

For Research Use Only

# SENP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83594-2-PBS

## Basic Information

<b>Catalog Number:</b> 83594-2-PBS	<b>GenBank Accession Number:</b> BC045639	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 29843	<b>CloneNo.:</b> 240478G6
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9POJ3	
<b>Isotype:</b> IgG	<b>Full Name:</b> SUMO1/sentrin specific peptidase 1	
<b>Immunogen Catalog Number:</b> AG18028	<b>Calculated MW:</b> 644 aa, 73 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00607-1, SENP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83594-2-PBS. Detection antibody: 83594-1-PBS. Standard: Ag18028. Range: 0.625-80 ng/mL