

For Research Use Only

# KMT2C Recombinant antibody, PBS Only (Detector)



Catalog Number:83595-2-PBS

## Basic Information

Catalog Number: 83595-2-PBS	GenBank Accession Number: NM_170606	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 58508	CloneNo.: 240542G9
Source: Rabbit	UNIPROT ID: Q8NEZ4	
Isotype: IgG	Full Name: myeloid/lymphoid or mixed-lineage leukemia 3	
Immunogen Catalog Number: AG28875	Calculated MW: 541 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00610-1, KMT2C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83595-1-PBS. Detection antibody: 83595-2-PBS. Standard: Ag28875. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00610-2, KMT2C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83595-3-PBS. Detection antibody: 83595-2-PBS. Standard: Ag28875. Range: 0.313-40 ng/mL