

For Research Use Only

BARHL2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83596-2-PBS

Basic Information

Catalog Number: 83596-2-PBS	GenBank Accession Number: BC126439	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 343472	CloneNo.: 240558F12
Source: Rabbit	UNIPROT ID: Q9NY43	
Isotype: IgG	Full Name: BarH-like homeobox 2	
Immunogen Catalog Number: AG21156	Calculated MW: 387 aa, 42 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00582-1, BARHL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83596-1-PBS. Detection antibody: 83596-2-PBS. Standard: Ag21156. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00582-3, BARHL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83596-3-PBS. Detection antibody: 83596-2-PBS. Standard: Ag21156. Range: 0.625-40 ng/mL