

For Research Use Only

BARHL2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83596-3-PBS

Basic Information

Catalog Number:

83596-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21156

GenBank Accession Number:

BC126439

GeneID (NCBI):

343472

UNIPROT ID:

Q9NY43

Full Name:

BarH-like homeobox 2

Calculated MW:

387 aa, 42 kDa

Purification Method:

Protein A purification

CloneNo.:

240558D3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00582-3, BARHL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83596-3-PBS. Detection antibody: 83596-2-PBS. Standard: Ag21156. Range: 0.625-40 ng/mL