

For Research Use Only

# Flotillin 1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83598-2-PBS

## Basic Information

Catalog Number: 83598-2-PBS	GenBank Accession Number: BC001146	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10211	CloneNo.: 240582E11
Source: Rabbit	UNIPROT ID: O75955	
Isotype: IgG	Full Name: flotillin 1	
Immunogen Catalog Number: AG7902	Calculated MW: 47 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00584-2, Flotillin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83598-3-PBS. Detection antibody: 83598-2-PBS. Standard: Ag7902. Range: 6.25-100 ng/mL.