

For Research Use Only

# CRY1 Recombinant antibody, PBS Only



Catalog Number: 83602-5-PBS

## Basic Information

Catalog Number: 83602-5-PBS	GenBank Accession Number: BC030519	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1407	CloneNo.: 240520B12
Source: Rabbit	UNIPROT ID: Q16526	
Isotype: IgG	Full Name: cryptochrome 1 (photolyase-like)	
Immunogen Catalog Number: AG4270	Calculated MW: 586 aa, 66 kDa	

## Applications

Tested Applications:  
WB, ELISA

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83602-5-RR (CRY1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83602-5-PBS in a different storage buffer formulation.