

For Research Use Only

# C11orf67 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83603-1-PBS

## Basic Information

<b>Catalog Number:</b> 83603-1-PBS	<b>GenBank Accession Number:</b> BC002752	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 28971	<b>CloneNo.:</b> 240530D4
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9H7C9	
<b>Isotype:</b> IgG	<b>Full Name:</b> chromosome 11 open reading frame 67	
<b>Immunogen Catalog Number:</b> AG14258	<b>Calculated MW:</b> 122 aa, 13 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00605-2, C11orf67 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83603-3-PBS. Detection antibody: 83603-1-PBS. Standard: Ag14258. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00605-3, C11orf67 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83603-1-PBS. Detection antibody: 83603-2-PBS. Standard: Ag14258. Range: 0.313-40 ng/mL.