

For Research Use Only

# C11orf67 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83603-3-PBS

## Basic Information

Catalog Number:

83603-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14258

GenBank Accession Number:

BC002752

GeneID (NCBI):

28971

UNIPROT ID:

Q9H7C9

Full Name:

chromosome 11 open reading frame

67

Calculated MW:

122 aa, 13 kDa

Purification Method:

Protein A purification

CloneNo.:

240530G5

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00605-1, C11orf67 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83603-3-PBS. Detection antibody: 83603-4-PBS. Standard: Ag14258. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00605-2, C11orf67 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83603-3-PBS. Detection antibody: 83603-1-PBS. Standard: Ag14258. Range: 0.313-40 ng/mL.