

For Research Use Only

C11orf67 Recombinant antibody, PBS Only (Detector)



Catalog Number:83603-4-PBS

Basic Information

Catalog Number:

83603-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14258

GenBank Accession Number:

BC002752

GeneID (NCBI):

28971

UNIPROT ID:

Q9H7C9

Full Name:

chromosome 11 open reading frame

67

Calculated MW:

122 aa, 13 kDa

Purification Method:

Protein A purification

CloneNo.:

240530H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, WB

Species Specificity:

rat, mouse, human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83603-4-RR (C11orf67 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83603-4-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP00605-1, C11orf67 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83603-3-PBS. Detection antibody: 83603-4-PBS. Standard: Ag14258. Range: 0.313-40 ng/mL.



Bi-layer interferometry (BLI) kinetic assays of 83603-4-RR against Human C11orf67 were performed. The affinity constant is 0.304 nM .