

For Research Use Only

POMP Recombinant antibody, PBS Only (Detector)

Catalog Number: 83607-2-PBS



Basic Information

Catalog Number: 83607-2-PBS	GenBank Accession Number: BC003390	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51371	CloneNo.: 240563D1
Source: Rabbit	UNIPROT ID: Q9Y244	
Isotype: IgG	Full Name: proteasome maturation protein	
Immunogen Catalog Number: AG6987	Calculated MW: 16 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00590-3, POMP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83607-3-PBS. Detection antibody: 83607-2-PBS. Standard: Ag6987. Range: 1.25-80 ng/mL