

For Research Use Only

ATM Recombinant antibody, PBS Only (Detector)



Catalog Number: 83608-3-PBS

Basic Information

Catalog Number: 83608-3-PBS	GenBank Accession Number: BC137169	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 472	CloneNo.: 240560H7
Source: Rabbit	UNIPROT ID: Q13315	
Isotype: IgG	Full Name: ataxia telangiectasia mutated	
Immunogen Catalog Number: AG25920		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00592-2, ATM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83608-1-PBS. Detection antibody: 83608-3-PBS. Standard: Ag25920. Range: 0.781-100 ng/mL.