

For Research Use Only

STEAP2 Recombinant antibody, PBS Only (Detector)

Catalog Number:83611-2-PBS



Basic Information

Catalog Number: 83611-2-PBS	GenBank Accession Number: NM_152999	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 261729	CloneNo.: 240611B11
Source: Rabbit	UNIPROT ID: Q8NFT2	
Isotype: IgG	Full Name: six transmembrane epithelial antigen of the prostate 2	
	Calculated MW: 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00600-1, STEAP2-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83611-4-PBS. Detection antibody: 83611-2-PBS. Standard: SY00725. Range: 0.039-10 ng/mL



Cytometric bead array standard curve of MP00600-2, STEAP2-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83611-1-PBS. Detection antibody: 83611-2-PBS. Standard: SY00725. Range: 0.039-10 ng/mL