

For Research Use Only

# STEAP2 Recombinant antibody, PBS Only (Capture)

Catalog Number:83611-4-PBS



## Basic Information

Catalog Number:

83611-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_152999

GeneID (NCBI):

261729

UNIPROT ID:

Q8NFT2

Full Name:

six transmembrane epithelial antigen of the prostate 2

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

240611F4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00600-1, STEAP2-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83611-4-PBS. Detection antibody: 83611-2-PBS. Standard: SY00725. Range: 0.039-10 ng/mL.