

For Research Use Only

STEAP2 Recombinant antibody, PBS Only (Detector)

Catalog Number:83611-5-PBS



Basic Information

Catalog Number:

83611-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_152999

GeneID (NCBI):

261729

UNIPROT ID:

Q8NFT2

Full Name:

six transmembrane epithelial antigen of the prostate 2

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

240611H3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00600-3, STEAP2-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83611-3-PBS. Detection antibody: 83611-5-PBS. Standard: SY00725. Range: 0.039-10 ng/mL.