

For Research Use Only

PTEN Recombinant antibody, PBS Only (Detector)



Catalog Number:83612-4-PBS

Basic Information

Catalog Number: 83612-4-PBS	GenBank Accession Number: BC005821	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5728	CloneNo.: 240521G12
Source: Rabbit	UNIPROT ID: P60484	
Isotype: IgG	Full Name: phosphatase and tensin homolog	
	Calculated MW: 47 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00328-4, phospho-PTEN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83612-1-PBS. Detection antibody: 83612-4-PBS. Standard: SY01845P. Range: 0.195-100 ng/mL.