

For Research Use Only

AP1M1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83613-3-PBS

Basic Information

Catalog Number: 83613-3-PBS	GenBank Accession Number: BC017469	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8907	CloneNo.: 240553A12
Source: Rabbit	UNIPROT ID: Q9BXS5	
Isotype: IgG	Full Name: adaptor-related protein complex 1, mu 1 subunit	
Immunogen Catalog Number: AG2757	Calculated MW: 423 aa, 49 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00588-3, AP1M1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83613-1-PBS. Detection antibody: 83613-3-PBS. Standard: Ag2757. Range: 3.125-100 ng/mL.