

For Research Use Only

AP1M1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83613-5-PBS

Basic Information

Catalog Number:

83613-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2757

GenBank Accession Number:

BC017469

GeneID (NCBI):

8907

UNIPROT ID:

Q9BXS5

Full Name:

adaptor-related protein complex 1,
mu 1 subunit

Calculated MW:

423 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

240553G9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00588-1, AP1M1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83613-2-PBS. Detection antibody: 83613-5-PBS. Standard: Ag2757. Range: 0.781-100 ng/mL.