

For Research Use Only

CNOT4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83615-2-PBS

Basic Information

Catalog Number: 83615-2-PBS	GenBank Accession Number: BC035590	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4850	CloneNo.: 240559C1
Source: Rabbit	UNIPROT ID: O95628	
Isotype: IgG	Full Name: CCR4-NOT transcription complex, subunit 4	
Immunogen Catalog Number: AG3163	Calculated MW: 575 aa, 64 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00593-3, CNOT4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83615-3-PBS. Detection antibody: 83615-2-PBS. Standard: Ag3163. Range: 0.313-40 ng/mL.