

For Research Use Only

Neurofibromin Recombinant antibody, PBS Only (Detector)



Catalog Number: 83620-3-PBS

Basic Information

Catalog Number: 83620-3-PBS	GenBank Accession Number: M60915	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4763	CloneNo.: 240540E12
Source: Rabbit	UNIPROT ID: P21359	
Isotype: IgG	Full Name: neurofibromin 1	
Immunogen Catalog Number: AG25799	Calculated MW: 319 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00570-3, Neurofibromin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83620-2-PBS. Detection antibody: 83620-3-PBS. Standard: Ag25799. Range: 0.313-40 ng/mL