

For Research Use Only

CDK8 Recombinant antibody, PBS Only (Detector)



Catalog Number:83621-4-PBS

Basic Information

Catalog Number: 83621-4-PBS	GenBank Accession Number: BC069634	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1024	CloneNo.: 240598C4
Source: Rabbit	UNIPROT ID: P49336	
Isotype: IgG	Full Name: cyclin-dependent kinase 8	
Immunogen Catalog Number: AG17148	Calculated MW: 464 aa, 53 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00571-2, CDK8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83621-3-PBS. Detection antibody: 83621-4-PBS. Standard: Ag17148. Range: 0.156-20 ng/mL