

For Research Use Only

MTMR15 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:83622-4-PBS



Basic Information

Catalog Number: 83622-4-PBS	GenBank Accession Number: BC047882	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 22909	CloneNo.: 240629D12
Source: Rabbit	UNIPROT ID: Q9Y2M0	
Isotype: IgG	Full Name: myotubularin related protein 15	
Immunogen Catalog Number: AG11651	Calculated MW: 60 kDa, 114 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00572-1, MTMR15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83622-4-PBS. Detection antibody: 83622-1-PBS. Standard: Ag11651. Range: 0.78-100 ng/mL.



Cytometric bead array standard curve of MP00572-3, MTMR15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83622-2-PBS. Detection antibody: 83622-4-PBS. Standard: Ag11651. Range: 0.78-100 ng/mL.