

For Research Use Only

# C1GALT1C1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83625-3-PBS

## Basic Information

Catalog Number: 83625-3-PBS	GenBank Accession Number: BC050441	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 29071	CloneNo.: 240624A10
Source: Rabbit	UNIPROT ID: Q96EU7	
Isotype: IgG	Full Name: C1GALT1-specific chaperone 1	
Immunogen Catalog Number: AG6899	Calculated MW: 36 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00609-1, C1GALT1C1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83625-4-PBS. Detection antibody: 83625-3-PBS. Standard: Ag6899. Range: 3.13-100 ng/mL.