

For Research Use Only

C1orf83 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83626-2-PBS

Basic Information

Catalog Number: 83626-2-PBS	GenBank Accession Number: BC035374	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 127428	CloneNo.: 240678B9
Source: Rabbit	UNIPROT ID: Q96MN5	
Isotype: IgG	Full Name: chromosome 1 open reading frame 83	
Immunogen Catalog Number: AG15040	Calculated MW: 208 aa, 23 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00578-1, C1orf83 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83626-1-PBS. Detection antibody: 83626-2-PBS. Standard: Ag15040. Range: 0.78-100 ng/mL.