

For Research Use Only

PD-1/CD279 Recombinant antibody, PBS Only (Capture)



Catalog Number:83632-2-PBS

Basic Information

Catalog Number:

83632-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074740

GeneID (NCBI):

5133

UNIPROT ID:

Q15116

Full Name:

programmed cell death 1

Calculated MW:

288 aa, 32 kDa

Purification Method:

Protein A purification

CloneNo.:

240724G7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00612-1, PD-1/CD279 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83632-2-PBS. Detection antibody: 83632-1-PBS. Standard: Eg0974. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00612-2, PD-1/CD279 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83632-2-PBS. Detection antibody: 83632-4-PBS. Standard: Eg0974. Range: 0.313-40 ng/mL.