

For Research Use Only

CDK2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83635-3-PBS

Basic Information

Catalog Number: 83635-3-PBS	GenBank Accession Number: BC003065	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1017	CloneNo.: 240580B9
Source: Rabbit	UNIPROT ID: P24941	
Isotype: IgG	Full Name: cyclin-dependent kinase 2	
Immunogen Catalog Number: AG17133	Calculated MW: 298 aa, 34 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00611-1, CDK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83635-3-PBS. Detection antibody: 83635-2-PBS. Standard: Ag17133. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00611-3, CDK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83635-3-PBS. Detection antibody: 83635-1-PBS. Standard: Ag17133. Range: 0.156-20 ng/mL