

For Research Use Only

ATP6V1E1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83636-2-PBS

Basic Information

Catalog Number: 83636-2-PBS	GenBank Accession Number: BC004443	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 529	CloneNo.: 240584H8
Source: Rabbit	UNIPROT ID: P36543	
Isotype: IgG	Full Name: ATPase, H ⁺ transporting, lysosomal 31kDa, V1 subunit E1	
Immunogen Catalog Number: AG7399	Calculated MW: 26 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00628-2, ATP6V1E1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83636-1-PBS. Detection antibody: 83636-2-PBS. Standard: Ag7399. Range: 0.625-40 ng/mL.