

For Research Use Only

SYT1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83637-2-PBS

Basic Information

Catalog Number: 83637-2-PBS	GenBank Accession Number: BC058917	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6857	CloneNo.: 240589H4
Source: Rabbit	UNIPROT ID: P21579	
Isotype: IgG	Full Name: synaptotagmin I	
Immunogen Catalog Number: AG5970	Calculated MW: 48 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00633-1, SYT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83637-2-PBS. Detection antibody: 83637-1-PBS. Standard: Ag5970. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00633-2, SYT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83637-2-PBS. Detection antibody: 83637-3-PBS. Standard: Ag5970. Range: 0.625-80 ng/mL