

For Research Use Only

# IFT20 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83639-3-PBS

## Basic Information

Catalog Number:

83639-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4521

GenBank Accession Number:

BC038094

GeneID (NCBI):

90410

UNIPROT ID:

Q8IY31

Full Name:

intraflagellar transport 20 homolog  
(Chlamydomonas)

Calculated MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

240552F3

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00644-1, IFT20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83639-3-PBS. Detection antibody: 83639-5-PBS. Standard: Ag4521. Range: 0.313-40 ng/mL