

For Research Use Only

IFT20 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83639-5-PBS

Basic Information

Catalog Number: 83639-5-PBS	GenBank Accession Number: BC038094	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 90410	CloneNo.: 240552F8
Source: Rabbit	UNIPROT ID: Q8IY31	
Isotype: IgG	Full Name: intraflagellar transport 20 homolog (Chlamydomonas)	
Immunogen Catalog Number: AG4521	Calculated MW: 15 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00644-1, IFT20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83639-3-PBS. Detection antibody: 83639-5-PBS. Standard: Ag4521. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00644-2, IFT20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83639-2-PBS. Detection antibody: 83639-5-PBS. Standard: Ag4521. Range: 0.313-40 ng/mL