

For Research Use Only

# KCND3 Recombinant antibody, PBS Only (Capture)

Catalog Number:83640-1-PBS



## Basic Information

Catalog Number:

83640-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22183

GenBank Accession Number:

BC113477

GeneID (NCBI):

3752

UNIPROT ID:

Q9UK17

Full Name:

potassium voltage-gated channel,  
Shal-related subfamily, member 3

Calculated MW:

655 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

240578A10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00616-2, KCND3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83640-1-PBS. Detection antibody: 83640-5-PBS. Standard: Ag22183. Range: 0.78-100 ng/mL